

PerfusionPal™ PP-A012-0001 PACKAGE CONTENTS

Item	Qty (Size)	Description
PerfusionPal™ insert system	1	Lid, insert, and tray with short, flexible tubing and a male luer fitting at the end of tubing
SeedEZ scaffolds	15 (0.5" dia.)	Glass fiber scaffolds for seeding cells
Seating insert	1	Specialized lid to be used to position SeedEZ scaffolds at the proper location in the wells.
Tubing	1 (6 ft.)	Silicone tubing with female luer fittings on both ends and a moveable stainless steel sheath for preventing crushing of tubing between incubator door and gasket.
Syringe	1 (60 mL)	Luer-lock fitting
Luer cap	1	Female luer cap to attach to male luer fitting on tray. <u>Purpose:</u> The luer cap plugs perfusion port for PerfusionPal use as an unperfused system.